



Series 204

Horizontal Laminar Flow Bench

The Series 204 laminar flow bench is a bench unit designed to install on top of typical metro wire or eagle-type wire shelving. This horizontal laminar flow bench unit is also designed to set on table top and blow horizontally.

Laminar flow bench units come in 2 x 2 ft., 2 x 3 ft., 2 x 4 ft., and 2 x 6 ft. sizes. The prefilter is serviced from behind the unit.

Features

- ❖ Decorative white painted steel housing
- ❖ Internal motorized filter unit for full filter coverage
- ❖ Clean hood on three sides
- ❖ Power cord

Options

- Stainless steel housing
- Stainless steel hood
- Light



For more information or to download or fax this product from the web, simply go to:

www.cleanairproducts.com/204

Specifications subject to change. Please contact factory for details.