



Series 416

Vertical Laminar Flow Clean Bench

The Series 416 is a vertical laminar flow clean bench unit designed to install on top of typical metro wire or eagle-type wire shelving.

The Series 416 is a vertical laminar flow bench unit designed to set on table top and blow vertically.

Units come in 2 x 2 ft., 2 x 3 ft., 2 x 4 ft., and 2 x 6 ft. sizes.

Features

- ❖ Decorative white painted steel housing
- ❖ Internal motorized filter unit for full filter coverage
- ❖ Solid curtains on three sides
- ❖ Strip door on front
- ❖ 2 ft. x 4ft. shelving size
- ❖ 120 V power cord

Series 416 Vertical Laminar Flow Clean Bench

Options

- Wire shelving
- Stainless steel housing
- Light
- Strip Door
- Hinged front doors
- Stainless steel sides

Guarantee

A written 1 year warranty is furnished with each clean bench.




Table Top Model

- Shown in 2x4 size
- White upper cabinet
- Flow thru light
- Stainless steel sides

Cart Model

- Shown in 2x4 size
- Stainless steel upper housing
- Curtains (3) sides
- Strip door front
- Stainless steel caster cart with solid bottom shelf



A = 24 in.
36 in.
48 in.
60 in.
72 in.

B* = 33.5 in.
53 in.
74 in.
86 in.
96 in.

C = 14 in.
20 in.**

* Post can be cut to custom height
**With flow thru light

Specifications subject to change. Please contact factory for details.

Clean Air
PRODUCTS

Solutions Built to Your Specifications.

8605 Wyoming Ave. N. • Minneapolis MN 55445

tel: 763.425.9122 800.423.9728 fax: 763.425.2004

e-mail: sales@cleanairproducts.com www.cleanairproducts.com