



NOTES: UNLESS OTHERWISE SPECIFIED

1. DIMENSIONS ARE IN INCHES (mm).




© COPYRIGHT 2010 CLEAN AIR PRODUCTS	Clean Air PRODUCTS	8605 WYOMING AVENUE NORTH • MINNEAPOLIS, MN 55445 800-423-9728 • ph: 763-425-9122 • fax: 763-425-2004 www.cleanairproducts.com			
	DESCRIPTION.	Horizontal Flow Clean Bench CAP301-830-30H			
	QUOTE NO.				
	DO NOT SCALE	DATE DRAWN 11/20/2012	QUOTE FORM 037723	REV 2	SHEET 1 of 2

NOTES: UNLESS OTHERWISE SPECIFIED

1. DIMENSIONS ARE IN INCHES (mm).

Options:



Standard Features:

1. Heavy Gauge White Vinyl Clad Steel Construction
2. White Plastic Laminate Table Top Material
3. Clear Polycarbonate Hood Panels
4. Fluorescent Light - T8 Lamps
5. Front Removable HEPA- 99.99% on 0.3 Micron Filter (Expanded Steel, Epoxy Painted, Screen On Final Filter Face)
6. Fiberglass Pre-filter Behind Hinged Louvered Grill
7. Electrical Service: Convenience Receptacle (4 Amps Available), 120V, 1-Phase, 3 Wire, With Power Cord, 17.38 A FLA
8. Blower & Light Switch
9. Power Cord (5-30P)
10. Table Top Height: 30" from Floor +/- 1/2" With Leg Adjustment
11. Blower Motor Assembly With Flex Duct
12. Adjustable Leg Leveling Glides

Options:

13. Magnehelic Gage
14. Minihelic Gage
15. Stainless Steel Work Surface
16. Stainless Steel Cabinet
17. Casters, With Equipment Legs - Increase Work Surface Height 6"
18. 6" Equipment Legs
19. Conductive Plastic Laminate Table Top
20. Ion Bar System
21. Grill on HEPA (Removable)
22. Anti-Static PVC Hood Panels
23. Safety Glass Hood Panels
24. Safety Glass Top & Stainless Hood Panels
25. Deeper Work Surface (27" & 33")
26. ULPA - 99.999% on 0.12 Micron Filter
27. Extra Service Receptacle-Separate 15A Power Cord
28. Remove Convenience Receptacle To Allow For A 5-20 Plug, This Would Achieve 13.38 A FLA

© COPYRIGHT 2010 CLEAN AIR PRODUCTS	Clean Air PRODUCTS 8605 WYOMING AVENUE NORTH • MINNEAPOLIS, MN 55445 800-423-9728 • ph: 763-425-9122 • fax: 763-425-2004 www.cleanairproducts.com	
	DESCRIPTION.	Horizontal Flow Clean Bench CAP301-830-30H
	QUOTE NO.	
	DO NOT SCALE	DATE DRAWN 11/20/2012 QUOTE FORM 037723 REV 2 SHEET 2 of 2