


1. DIMENSIONS ARE IN INCHES (mm).




Clean Air PRODUCTS	8605 WYOMING AVENUE NORTH • MINNEAPOLIS, MN 55445 800-423-9728 • ph: 763-425-9122 • fax: 763-425-2004 www.cleanairproducts.com				
DESCRIPTION.	Horizontal Flow Clean Bench CAP301-736-36H				
QUOTE NO.	_____				
DO NOT SCALE	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 33%; text-align: center;">DATE DRAWN 03/12/2021</td> <td style="width: 33%; text-align: center;">QUOTE FORM 115879</td> <td style="width: 16%; text-align: center;">REV 1</td> <td style="width: 18%; text-align: center;">SHEET 1 of 2</td> </tr> </table>	DATE DRAWN 03/12/2021	QUOTE FORM 115879	REV 1	SHEET 1 of 2
DATE DRAWN 03/12/2021	QUOTE FORM 115879	REV 1	SHEET 1 of 2		

NOTES: UNLESS OTHERWISE SPECIFIED

1. DIMENSIONS ARE IN INCHES (mm).

Options:



Standard Features:

1. Heavy Gauge White Vinyl Clad Steel Construction
2. White Plastic Laminate Table Top Material
3. Clear Polycarbonate Hood Panels
4. Fluorescent Light - T8 Lamps
5. Front Removable HEPA- 99.99% on 0.3 Micron Filter (Expanded Steel, Epoxy Painted, Screen On Final Filter Face)
6. Fiberglass Pre-filter Behind Hinged Louvered Grill
7. Electrical Service: Convenience Receptacle (4 Amps Available), 120V, 1-Phase, 3 Wire, With Power Cord, 16.89 A FLA
8. Blower & Light Switch
9. Power Cord (5-30P)
10. Table Top Height: 30" from Floor +/- 1/2" With Leg Adjustment
11. Blower Motor Assembly With Flex Duct
12. 6" Equipment Legs

Options:

13. Magnehelic Gage
14. Stainless Steel Work Surface
15. Casters, With Equipment Legs
16. Conductive Plastic Laminate Table Top
17. Ion Bar System
18. Grill on HEPA (Removable)
19. Anti-Static Polycarbonate Hood Panels
20. Safety Glass Hood Panels
21. Safety Glass Top & Stainless Hood Panels
22. Deeper Work Surface (27-1/2" & 33-1/2")
23. ULPA - 99.999% on 0.12 Micron Filter
24. Extra Service Receptacle-Separate 15A Power Cord

© COPYRIGHT 2010 CLEAN AIR PRODUCTS	Clean Air PRODUCTS	8605 WYOMING AVENUE NORTH • MINNEAPOLIS, MN 55445 800-423-9728 • ph: 763-425-9122 • fax: 763-425-2004 www.cleanairproducts.com			
	DESCRIPTION.	Horizontal Flow Clean Bench CAP301-736-36H			
	QUOTE NO.				
	DO NOT SCALE	DATE DRAWN 03/12/2021	QUOTE FORM 115879	REV 1	SHEET 2 of 2