




NOTES: UNLESS OTHERWISE SPECIFIED

1. DIMENSIONS ARE IN INCHES (mm).



|                                        |                                                                                                                                                                |                          |                      |                 |
|----------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------|----------------------|-----------------|
| © COPYRIGHT 2010<br>CLEAN AIR PRODUCTS | <b>Clean Air PRODUCTS</b> 8605 WYOMING AVENUE NORTH • MINNEAPOLIS, MN 55445<br>800-423-9728 • ph: 763-425-9122 • fax: 763-425-2004<br>www.cleanairproducts.com |                          |                      |                 |
|                                        | DESCRIPTION. Horizontal Flow Clean Bench<br>CAP303-636-36H                                                                                                     |                          |                      |                 |
|                                        | QUOTE NO.                                                                                                                                                      |                          |                      |                 |
|                                        | <b>DO NOT SCALE</b>                                                                                                                                            | DATE DRAWN<br>04/15/2021 | QUOTE FORM<br>116830 | REV<br>1        |
|                                        |                                                                                                                                                                |                          |                      | SHEET<br>1 of 2 |

## OPTIONS:



### Standard Features:

1. Heavy Gauge White Vinyl Clad Steel Construction
2. Clear Polycarbonate Hood Panels: 30" from Floor +/- 1/2" With Leg Adjustment
3. Front Removable HEPA- 99.99% on 0.3 Micron Filter (Expanded Steel, Epoxy Painted, Screen On Final Filter Face)
4. Fiberglass Pre-filter Behind Hinged Louvered Grill
5. Electrical Service: Convenience Receptacle (4 Amps Available), 120V, 1-Phase, 3 Wire, With Power Cord, 16.89 A FLA
6. Blower & Light Switch
7. Power Cord (5-30P)
8. Blower Motor Assembly With Flex Duct
9. 6" Equipment Legs

### Options:

10. Fluorescent Light - T8 Lamps
11. Magnehelic Gage
12. Casters, With Equipment Legs
13. Conductive Plastic Laminate Table Top
14. Ion Bar System
15. Grill on HEPA (Removable)
16. Anti-Static Polycarbonate Hood Panels
17. Safety Glass Hood Panels
18. Safety Glass Top & Stainless Hood Panels
19. Deeper Hood Area (36", 42", 48")
20. ULPA - 99.999% on 0.12 Micron Filter
21. Extra Service Receptacle-Separate 15A Power Cord

NOTES: UNLESS OTHERWISE SPECIFIED

1. DIMENSIONS ARE IN INCHES (mm).

|                                        |                                                                                                                                                                |                          |                      |                             |
|----------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------|----------------------|-----------------------------|
| © COPYRIGHT 2010<br>CLEAN AIR PRODUCTS | <b>Clean Air PRODUCTS</b> 8605 WYOMING AVENUE NORTH • MINNEAPOLIS, MN 55445<br>800-423-9728 • ph: 763-425-9122 • fax: 763-425-2004<br>www.cleanairproducts.com |                          |                      |                             |
|                                        | DESCRIPTION. Horizontal Flow Clean Bench<br>CAP303-636-36H                                                                                                     |                          |                      |                             |
|                                        | QUOTE NO.                                                                                                                                                      |                          |                      |                             |
|                                        | <b>DO NOT SCALE</b>                                                                                                                                            | DATE DRAWN<br>04/15/2021 | QUOTE FORM<br>116830 | REV<br>1<br>SHEET<br>2 of 2 |