


NOTES: UNLESS OTHERWISE SPECIFIED  
1. DIMENSIONS ARE IN INCHES (mm).



**Features:**

1. Clear Tempered Safety Glass
2. Stainless Steel Lift Off Hinges
3. Stainless Steel Door Frame
4. Over Center Compression Latches With Linkage Bar - Stainless Steel

**Options Shown:**


5. One Step, Pin Style, Mechanical Interlock
6. 2 Pc. Flush Mount Stainless Steel Mounting Frame Set

**Notes:**

- Cabinet Material: Stainless Steel, Type 304, #4 Finish
- Electrical Requirements: None
- Interlocking and Non-interlocking Units Have The Same Exterior Dimensions Unless Noted.
- Wall Cutout is  **$24\text{-}13/16"$  x  $23\text{-}11/16"$**   
( $[A] + 1/2"$  x  $[B] + 1/2"$ )

**Pass Thru Dimensions:**

**W = Inside Width**  
**H = Inside Height**  
**D = Inside Depth**  
**A = W + 6 5/16"**  
**B = H + 5 3/16"**



|                                        |                                                                                                                                                                |                         |                      |          |                 |
|----------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------|----------------------|----------|-----------------|
| © COPYRIGHT 2019<br>CLEAN AIR PRODUCTS | <b>Clean AiR PRODUCTS</b> 8605 WYOMING AVENUE NORTH • MINNEAPOLIS, MN 55445<br>800-423-9728 • ph: 763-425-9122 • fax: 763-425-2004<br>www.cleanairproducts.com |                         |                      |          |                 |
|                                        | DESCRIPTION. <b>Pass Thru - Wall Mount, Bio-Design, Stainless Steel - CAP18WBD-SST-18Wx18Hx18D</b>                                                             |                         |                      |          |                 |
|                                        | QUOTE NO.                                                                                                                                                      |                         |                      |          |                 |
|                                        | <b>DO NOT SCALE</b>                                                                                                                                            | DATE DRAWN<br>4/16/2020 | QUOTE FORM<br>108567 | REV<br>1 | SHEET<br>1 of 1 |